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~Ceramics Cluster Filler for Innovation~

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Coronas provisionales, etc.

Índice

- 3 Ceramics Cluster Technology
- 4 Test Result Data for Physical Properties
- 5 Reproducing Natural Tooth Appearance
- 6 TWiNY Flow
- 7-8 TWiNY Enamel Clear
- 9 TWiNY Basic Layering Steps
- 10 Basic Build-up Layering for Jacket Crown
- 11 Basic Build-up Layering for Inlay and Onlay
- 12 Polishing TWiNY
- 13-14 Shade Color Table
- 15-16 TWiNY Product Line-up
- 17 Yamakin Brand

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
Steps for Polishing TWiNY and TWiNY Flow

Diamond Polisher	Abrasive Particle Size of diamond (μm)	Rotation Speed Rate (rpm)
C&B Diamond Polisher	1 ~ 2	Lower than 10,000 rpm.
C&B Nano Diamond Polisher	0.5 or smaller	Lower than 20,000 rpm.

Related Products

C&B Diamond Polisher
8 g

C&B Nano Diamond
5 g

Shade Color Table

Basic Shades															
Invisible Opaque 2ml															IvO
Opaque 2ml															OA1 OA2 OA3 OA3.5 OA4 OB1 OB2 OB3 OB4 OC1 OC2 OC3 OC4 OD2 OD3 OD4
Cervical 2.6ml															CA1 CA2 CB1 CB2 CC1 CC2 CD1 CD2
Opaque Dentine 2.6ml															ODA1 ODA2 ODA3 ODA3.5 ODA4 ODB1 ODB2 ODB3 ODB4 ODC1 ODC2 ODC3 ODC4 ODD2 ODD3 ODD4
Dentine 2.6ml															DA1 DA2 DA3 DA3.5 DA4 DB1 DB2 DB3 DB4 DC1 DC2 DC3 DC4 DD2 DD3 DD4
Enamel 2.6ml															E2 E3 E4 E1 E2 E3 E2 E3 E4 E3 E4
Enamel Clear															
4.8g(2.6ml)															
Gum Shades															
Gum Opaque 2ml															
Gum 2.6ml															
Gum Modifier 2.6ml															
Opaque Special Colors 2ml															
Effect 2.6ml															
Translucent 2.6ml															
TWiNY Flow 2ml															
Note: DA1 Flow, DA2 Flow, DA3.5 Flow, DA4 Flow, DB1 Flow, E1 Flow, E2 Flow, E3 Flow, CA1 Flow, CA2 Flow															
DA1 Flow DA2 Flow DA3 Flow DA3.5 Flow DA4 Flow DB1 Flow E1 Flow E2 Flow E3 Flow CA1 Flow CA2 Flow															
Coffee Flow Orange Flow AM Flow OC Flow WE Flow															
TE Flow HVT Flow LVT Flow CT1 Flow CT2 Flow CT3 Flow CT4 Flow T Glass Flow T Glass Clear Flow T Blue Flow															
G1 Flow G3 Flow G5 Flow G7 Flow G Or Flow GM Gray Flow GM Trans Flow															
Others															
7ml Resin Air Barrier															
5ml Resin Separator															
5ml Resin Spacer															

Stain

White	Milky	Pink	Salmon Pink	White Violet	Violet	Orange	Dark Orange	Red
Yellow	Blue	Brown	Reb Brown	Dark Brown	A Shift	B Shift	C Shift	D Shift
Gray	Black	Clear						

Luna-Wing stains can be used for TWiNY (internal use).

Basic Shades

Invisible Opaque
Invisible Opaque is low-flowable opaque resin that flows into the undercut of retention beads to enhance physical bonding strength between metal and resin.

Opaque
Opaque resin for metal surface control.

Opaque Dentine
Body resin to prevent percolation of opaque color when there is extremely thin space for build-up, and to express further color depth. Apply thinly under Dentine layer when there is 0.5 mm thickness.

Cervical
Body resin for natural color expression around cervical area.

Dentine
Body resin to express dentine.

Enamel
Body resin to express enamel.

[Enamel Clear]

Enamel to express higher translucency than usual Enamel.

[Trans Enamel]

To express transparency or translucency between Translucent and Enamel.

[Translucent]

Body resin to express translucency.

HVT	High Value Translulent	Translucent with high brightness
T	Translucent	Normal translucent
LVT	Low Value Translulent	Translucent with low brightness
CT1		Light orange and pink translucent for gum color. For cervical area expression.
CT2	Cervical Translulent	Light orange and pinkish translucent for cervical area.
CT3		Orange translucent for cervical area.
CT4		Yellowish translucent for cervical area.
T Glass	Translucent Glass	The second highest transparency of TWiNY Line-up.
T Glass Clear	Translucent Glass Clear	Translucent with highest transparency in the lineup.
T Blue	Translucent Blue	Blue Translucent for incisal edge.

[Effect]

Characterizing Body resin to express discolored teeth.

HV WE	High Value White Enamel	Enamel Color with higher intensity (whiter) than WE.
WE	White Enamel	Enamel color with high intensity. To express proximal surface and white bands. Whiter than EO and not as transparent as Enamel EO.
AM	Amber	To express orangish translucent (amber color).
OC	Occlusal	Very light orange color for occlusal area of molars.
Coffee	Coffee	To accent brown color like coffee.
Orange	Orange	To accent orange color.

[Base]

Body resin as base to fill pontic part of bridge, etc.

Translucent color with deep light-cure depth designed only for pontic part. It cannot be used on facing part.

[TWiNY Flow]

TWiNY Flow is useful for cases such as applying to cervical area and cavity bottoms, the color expression of pits and fissures, building up Enamel ridges, adjusting color tones of gingival area, repairing air bubbles, pouring in pontic area, etc.

[Gum Shades]

etc.

OG1	Gum Opaque	To express light pink gum.
OG2		To express orange pink gum.
OG3		To express normal orange pink gum.
OG4		To express dark orange pink gum.
OG5		To express dark red gum.
OG Or		To express orange gum.
G1	Gum	Light pink gum color to express gum.
G2		Dark pink gum color to express gum.
G3		Light orange pink gum color to express gum.
G4		Normal gum color to express gum.
G5		Orange pink gum color to express gum.
G6		Dark orange pink gum color to express gum.
G7		Dark red gum color to express gum.
G Or	Gum	Orange gum color to express gum.
GM Gray		To express melanin pigment.
GM Trans	Gum Modifier	To express transparency on the surface of gum.

[Gum Stain]

Characterizing paste to express discolored teeth. Expresses effective color tone by thin layer application. Apply with a round brush after the light

TWiNY Product Line-up

ADVANCED SET

Invisible Opaque 2mL : IvO
Opaque 2mL : OA2, OA3, OA3.5
Opaque Special Colors 2mL : InO1, MO
Cervical 4.8g(2.6mL) : CA1, CA2
Opaque Dentine 4.8g(2.6mL) : ODA2, ODA3, ODA3.5
Accessory • Flat Brush 3 pcs • Round Brush 3 pcs • Mixing Papers 50 sheets • Palette 5 pcs (1 shade cover)

Dentine 4.8g(2.6mL) : DA2, DA3, DA3.5
Enamel 4.8g(2.6mL) : E3, E4
Translucent 4.8g(2.6mL) : T, CT2
Gum Opaque 2mL : OG1, OG3, OG5
Gum Modifier 4.8g(2.6mL) : GM Gray, GM Trans
Gum Stain 1mL : G Dark Red, G Red, G Milky, G Violet
Multi Primer Paste 2mL
Multi Primer Repair Liquid One 6mL
Resin Spacer 5mL


REGULAR SET

Invisible Opaque 2mL : IvO
Opaque 2mL : OA2, OA3, OA3.5, OB2, OB3
Opaque Special Colors 2mL : InO1, MO
Cervical 4.8g(2.6mL) : CA1, CA2, CB1
Opaque Dentine 4.8g(2.6mL) : ODA2, ODA3, ODA3.5, ODB2, ODB3
Accessory • Flat Brush 3 pcs • Round Brush 3 pcs • Mixing Papers 50 sheets • Palette 5 pcs (1 shade cover)

Dentine 4.8g(2.6mL) : DA2, DA3, DA3.5, DB2, DB3
Enamel 4.8g(2.6mL) : E2, E3, E4
Translucent 4.8g(2.6mL) : T, CT2
Multi Primer Paste 2mL
Multi Primer Repair Liquid One 6mL
Resin Spacer 5mL

INTRO SET

Invisible Opaque 2mL : IvO
Opaque 2mL : OA3
Cervical 4.8g(2.6mL) : CA1
Dentine 4.8g(2.6mL) : DA3
Enamel 4.8g(2.6mL) : E3
Translucent 4.8g(2.6mL) : T
Multi Primer Paste 2mL
Multi Primer Repair Liquid One 6mL
Accessory • Flat Brush 2 pcs • Round Brush 2 pcs

Multi Primer

KZR-CAD HR 2 GR

Hybrid Ceramics Block for CAD/CAM Use

KZR-CAD HR 2

Hybrid Ceramics Block for CAD/CAM Use

Luna-Wing

Light-cured Indirect Composite Resin

Multi Primer

Bonding Materials

STARTER SET

Invisible Opaque 2mL : IvO
Opaque 2mL : OA2, OA3, OA3.5
Opaque Special Colors 2mL : InO1, MO
Cervical 4.8g(2.6mL) : CA1, CA2
Opaque Dentine 4.8g(2.6mL) : ODA2, ODA3, ODA3.5
Accessory • Flat Brush 3 pcs • Round Brush 3 pcs • Mixing Papers 50 sheets • Palette 5 pcs (1 shade cover)

Opaque Dentine 4.8g(2.6mL) : ODA2, ODA3, ODA3.5
Enamel 4.8g(2.6mL) : E3, E4
Translucent 4.8g(2.6mL) : T, CT2
Multi Primer Paste 2mL
Multi Primer Repair Liquid One 6mL
Resin Spacer 5mL

GUM KIT

Gum Opaque 2mL : OG1, OG2, OG3, OG4, OG5
Gum 4.8g(2.6mL) : G1, G2, G3, G4, G5, G6, G7
Gum Modifier 4.8g(2.6mL) : GM Gray, GM Trans
Gum Stain 1mL : G Dark Red, G Red, G Milky, G Violet
Accessory • Flat Brush 3 pcs • Round Brush 3 pcs • Mixing Papers 50 sheets • Palette 5 pcs (1 shade cover)